

**Advanced Research Biological Microscope
CPD.43.SC53-3C026**



CPD.43.SC53-3C026				
MAIN BODY				
Item	Model	PN	Descriptions	QTY
10	SC53F	N002325	SC53 FRAME FORR/T. LIGHT	1
20	U-LHLEDC10	N002347	TRANSMITTED LAMP HOUSING	1
30	SC3-ARM	N002328	ARM FOR SC53F AND SC63F FRAME	1
TUBE&EYEPIECE				
40	U-TR30-2	N002330	TRINOCULAR TUBE, 0/100,50/50	1
50	SWH10X-H	N002253	WIDEFIELD EYEPIECE 10X, F.N=26.5FOCUSABLE	2
STAGE & CONDENSER				
60	U-SVRB	N002329	2-SLIDES MECHANICAL STAGE	1
70	U-HLDT	N001317	SPECIMEN HOLDER FORMECHANICAL STAGE	1
80	U-SC3	N000909	SWING-OUT CONDENSER N.A. 0.9	1
OBJECTIVE & NOSEPIECE				
90	U-D6RE	N001598	SEXTUPLE REVOLVING NOSEPIECE-RMS	1
100	UPLN4X	N002221	PLANACHROMAT OBJECTIVE 4X/0.13 WD 17	1
110	UPLN10X	N002222	PLAN ACHROMAT OBJECTIVE 10X/0.30 WD 10	1
120	UPLN20X	N002223	PLAN ACHROMAT OBJECTIVE 20X/0.50 WD 3.0	1
130	UPLN40X	N002225	PLAN ACHROMAT OBJECTIVE 40X/0.75 WD 0.6	1
140	UPLN100XO	N002229	PLAN ACHROMAT OBJECTIVE 100X/1.30 WD 0.2	Optional
ACCESSORIES				
150	45S-LBA	N001770	LBA FILTER	1
160	D630	N002022	6.3MP DIGITAL CAMERA USB3.0 WITH SOFTWARE	Optional
170	U-TV0.5XC	N001020	C-MOUNT VIDEO ATTACHMENT WITH 0.5X RELAY LENS	1
180	COVER-03	N000776	DUST COVER	1
190	UYCP-15	N001397	POWER CORD	1

